




CARFIL QUALITY
YOUR SUPPLIER FOR ALL LABORATORY NEEDS

CSF Sampling Port



Features

For CSF Sampling in large animals, use the MIN-CSF-SAM, a specialized port device which provides a subcutaneous guide cannula which is accessed via the port septum using a standard spinal needle.

The schematic above shows the device's cranial placement which is accomplished with stereotactic technique in order to orient the guide cannula to the proper CSF space.



PAVAN SERVICE BVBA

Beyntellus 3 - B-2360 Oud-Turnhout - Tel. +32 (0)14 45 13 10 - admin@carfil.be - BTW BE 0420.681.872

www.carfil.be