



CARFIL QUALITY
YOUR SUPPLIER FOR ALL LABORATORY NEEDS

Gas Anesthesia Mask



Anesthetics are required in conducting the experiment with the Stereotaxic Instrument, thus the anesthesia mask is used to match with the anesthesia machine, which prevents the researcher's exposure to the anesthetic gas. Together with the Gas Anesthesia Mask, the platform provides precise and secure positioning for mouse or rat anesthesia and surgery. This mask is installed in the adaptor, with the cone opening towards the "U" opening.

Order No.	Product information
68U601	Gas Anesthesia Mouse Mask
68U602	Gas Anesthesia Rat Mask



PAVAN SERVICE BVBA

Beyntellus 3 - B-2360 Oud-Turnhout - Tel. +32 (0)14 45 13 10 - admin@carfil.be - BTW BE 0420.681.872

www.carfil.be