



CARFIL QUALITY
YOUR SUPPLIER FOR ALL LABORATORY NEEDS

68U101 Microdriver (10 μ m)



In electrophysiological stimulation or recording experiments, it requires critical accurate electrode manipulation. This microdriver will improve operational accuracy in the vertical direction to 10 μ m. It is easy to handle and various electrode holders are available to work with the microdriver.



PAVAN SERVICE BVBA

Beyntellus 3 - B-2360 Oud-Turnhout - Tel. +32 (0)14 45 13 10 - admin@carfil.be - BTW BE 0420.681.872

www.carfil.be