



68U398 A/P Zero Point Calibrator



This A/P Zero appliance can be set electrodes position on the large animal Stereotaxic Instruments. It includes an A/P rail of 255mm long, with the scale of +100mm~ at 1m accuracy. This rail is fixed to an aluminum plate (255x128mm), and there is a small cyliner aligned to the zero point of the rail, which is coordinated with the zero point of the ear bars and the rail on the stereotaxic. Combined with the 2-D manipulator and the A/P rail, this zero appliance can calibrate the

electrode to Zero position in a three-dimensional coordinate system.



