



CARFIL QUALITY
YOUR SUPPLIER FOR ALL LABORATORY NEEDS

68U398 A/P Zero Point Calibrator



This A/P Zero appliance can be set electrodes position on the large animal Stereotaxic Instruments. It includes an A/P rail of 255mm long, with the scale of +100mm~ at 1mm accuracy. This rail is fixed to an aluminum plate (255x128mm), and there is a small cylinder aligned to the zero point of the rail, which is coordinated with the zero point of the ear bars and the rail on the stereotaxic. Combined with the 2-D manipulator and the A/P rail, this zero appliance can calibrate the electrode to Zero position in a three-dimensional coordinate system.



PAVAN SERVICE BVBA

Beyntellus 3 - B-2360 Oud-Turnhout - Tel. +32 (0)14 45 13 10 - admin@carfil.be - BTW BE 0420.681.872

www.carfil.be