



CARFIL QUALITY
YOUR SUPPLIER FOR ALL LABORATORY NEEDS

68U606 Manual Syringe Pump



The Manual Stereotaxic Injector provides simple, precise manual control to inject/withdraw the desired amount of fluid directly through the needle in the brain. It can be mounted directly to a stereotaxic manipulator and the easy control knob allows the fine movement of the micro-syringe. The total range of movement is 25mm with accuracy of 10 μ m. The syringe volume ranges from 1 μ l to 100 μ l.

(syringe not included)



PAVAN SERVICE BVBA

Beyntellus 3 - B-2360 Oud-Turnhout - Tel. +32 (0)14 45 13 10 - admin@carfil.be - BTW BE 0420.681.872

www.carfil.be