



CARFIL QUALITY
YOUR SUPPLIER FOR ALL LABORATORY NEEDS

78U001 Microdrill



The high speed drill is designed for opening the skull in stereotaxic surgery, to implant electrodes, canula, microdialysis probes, etc. This Drill is easy to handle and operate.

Technical Specifications

- Speed up to 38,000rpm
- Available options: forward or reverse rotation.
- Can be controlled manually or by foot paddle.
- Can be installed to the stereotaxic by the drill holder (68U605) and controlled by the manipulator.
- Various drills available, size from 0,5 -2,3mm in diameter.

Order No.	Product description
78U001	Micro Drill (Drill bits not included)
68U605	Microdrill holder
78U041	Drill bits, 0,6 mm HM1006 (set of 5pcs) Round tip
78U042	Drill bits, 0,8 mm HM1008 (set of 5pcs) Round tip
78U043	Drill bits, 1,0 mm HM1010 (set of 5pcs) Round tip
78U044	Drill bits, 1,2 mm HM1012 (set of 5pcs) Round tip
78U045	Drill bits, 1,4 mm HM1014 (set of 5pcs) Round tip
78U046	Drill bits, 1,6 mm HM1016 (set of 5pcs) Round tip
78U047	Drill bits, 1,8 mm HM1018 (set of 5pcs) Round tip
78U050	Drill bits, 2,3 mm HM1023 (set of 5pcs) Round tip
78U060	Drill bits, 2,7 mm HM1027 (set of 5pcs) Round tip
78U040	Drill bits, 0,5 mm HM1005 (set of 5pcs) Round tip



PAVAN SERVICE BVBA

Beyntellus 3 - B-2360 Oud-Turnhout - Tel. +32 (0)14 45 13 10 - admin@carfil.be - BTW BE 0420.681.872

www.carfil.be



CARFIL QUALITY
YOUR SUPPLIER FOR ALL LABORATORY NEEDS

78U061	Drill bits, 0,8 mm HM31008 (set of 5pcs) Flat tip
78U062	Drill bits, 1,0 mm HM31010 (set of 5pcs) Flat tip
78U063	Drill bits, 1,2 mm HM31012 (set of 5pcs) Flat tip
78U064	Drill bits, 1,4 mm HM31014 (set of 5pcs) Flat tip
78U065	Drill bits, 1,6 mm HM31016 (set of 5pcs) Flat tip
78U066	Drill bits, 1,8 mm HM31018 (set of 5pcs) Flat tip
78U067	Drill bits, 2,1 mm HM31021 (set of 5pcs) Flat tip
78U068	Drill bits, 2,3 mm HM31023 (set of 5pcs) Flat tip



PAVAN SERVICE BVBA

Beyntellus 3 - B-2360 Oud-Turnhout - Tel. +32 (0)14 45 13 10 - admin@carfil.be - BTW BE 0420.681.872

www.carfil.be