

4-way valves



Caractéristiques

To be used for simultaneously changing the direction of both the **SUPPLY** of the anaesthetic gas mixture **AND** the **ECHAUST** of the anaesthetic gas mixture. The risk of making mistakes when the gas is supply and gas exhaust have to be changed individually from f.i. the Induction box to the face mask, is reduced substantially bij using these valves.

4-Way Valve

A.. By turning the knob to the left, both supply and exhaust are being directed to/from the induction box.

B.. By turning the knob to the center, both supply and exhaust are being directed to/form the intubation aid.

C.. By turning the knob to the right, both supply and exhaust are being directed to/from the UMV-03

Numéro de commande

180000259

Double 4-way valve



