

Controlled Heating System Control Unit-MS/CU-MS



Control Unit-M/S

The **control Unit M/s** can like the Control Unit 02 be connected to a heating device and a temperature probe. However the interaction between the 2 connections allows you to **set the required animals rectal temperature**. The rectal probe measures the animals temperature and feed back is given to the controles of the heating device to keep the animals temperature as constant as possible. A safety feature is built-in to prevent over-heating in case the rectal probe is not in the animal. The heating element will not heat up over 42?C (also an acceptable norm for human skin).

With the Control Unit M/S you have a **self-regulating system** and do not need to manually adjust the temperature settings during the procedure.

Ordering number

Bestelnummer	
180000124	Control Unit MS
180000028	Flat Plate / 34,5 x 23 cm
180000128	Flat Plate / 18,5 x 11 cm
180000027	Half-Pipe Heating Unit
180000029	Tail Heating Unit
180000066	Temperature Probe - Mouse
180000026	Temperature Probe - Rat