



CARFIL QUALITY
YOUR SUPPLIER FOR ALL LABORATORY NEEDS

Pegasus pumps



	Pegasus 500C	Pegasus 500L
	Micropulsing Pump	Peristaltic
Flowrates	0,01 - 15,00 ml/h	0,10 - 100 ml/h
Pulses	0,8 µl	7 µl
Size	8,4*6,2*3,1 cm	8,4*6,2*3,1 cm
Weight	180 g	180 g
Accuracy	+/- 2%	+/- 5%

Features & Specifications

- Small Size
- Low weight
- Quiet Operation
- Battery Efficient
- High accurate flow

Drugbags

Drug bag with luer lock connector

- Filling using a syringe
- Available in different sizes - 50 ml, 100 ml, 150 ml and 300 ml
- Latex-free

Pegaboxes



PAVAN SERVICE BVBA

Beyntellus 3 - B-2360 Oud-Turnhout - Tel. +32 (0)14 45 13 10 - admin@carfil.be - BTW BE 0420.681.872

www.carfil.be



CARFIL QUALITY
YOUR SUPPLIER FOR ALL LABORATORY NEEDS

For holding and protection of the drugbag. Available in three different sizes (50 ml, 100 ml and 150 ml).

Pumpreader

The PEGA® PUMP READER SET comprises the PEGA® INTERFACE 2.0 for connecting the PEGA® infusion pump to a computer, and documentation software, the PEGA® PUMP READER.

Fully integratable in the Orchestra System.



PAVAN SERVICE BVBA

Beyntellus 3 - B-2360 Oud-Turnhout - Tel. +32 (0)14 45 13 10 - admin@carfil.be - BTW BE 0420.681.872

www.carfil.be