



CARFIL QUALITY
YOUR SUPPLIER FOR ALL LABORATORY NEEDS

68U062 Rat Gas Anesthesia Head Holder



This design integrates the gas anesthesia mask function while holding the rats head firmly. The adaptor can rotate the head to multi-angels without blocking the anesthetic gas delivery. The Dorsal/Ventral adjustment rage is 40° ($+20^{\circ} \sim -20^{\circ}$) with 1° accuracy. The horizontal adjustment range is 43,5mm, and the diameter of the inlet/outlet of gas anesthesia mask is 2mm.

Please make sure to use a safe gas evacuation!



PAVAN SERVICE BVBA

Beyntellus 3 - B-2360 Oud-Turnhout - Tel. +32 (0)14 45 13 10 - admin@carfil.be - BTW BE 0420.681.872

www.carfil.be