



## 68U066 Neonatal Rat Gas Anesthesia Head Holder



This model offers maximum stereotaxic adjustment for positioning the head, allowing versatility to angle the head without obstructing the delivery of the anesthetic gas. For different aged neonatal rat, there are 4 options for the gas mask: 10mm, 10,5mm, 11mm and 11,5mm in diameter.

- The vertical adjustable range: 30 mm (+ $10 \text{mm} \sim -20 \text{mm}$ ), with  $100 \mu \text{m}$  accuracy.
- The vertical rotatable range:  $70^{\circ}$  (-35mm ~ +35mm), with  $5^{\circ}$  accuracy.
- The horizontal rotatable range:  $70^{\circ}$  (-35mm  $\sim$  +35mm), with  $5^{\circ}$  accuracy.
- The horizontal adjustable range: 43,5mm.
- The diameter of the inlet/outlet ports for the mask is 2mm.

Please make sure to use a safe gas evacuation!