



68U061 Mouse Gas Anesthesia Head Holder



- This adaptor offers maximum stereotaxic adjustment for positioning the head allowing versatility to angle the head without obstructing the delivery of the anesthetic gas.
- The vertical adjustable range: 30mm (+10mm \sim -20mm), with 100 μ m accuracy.
- The vertical rotatable range: 70° (-35mm $\sim +35$ mm), with 5° accuracy.
- The horizontal rotatable range: 70° (-35mm $\sim +35$ mm), with 5° accuracy.
- The horizontal adjustable range: 43,5mm.
- The diameter of the inlet/outlet ports for the mask: 2m