



CARFIL QUALITY
YOUR SUPPLIER FOR ALL LABORATORY NEEDS

Canulla Holders



68U205

This holder is suitable for 3,5mm cannula heads as well as that of the micro-syringe pump.

68U203

The wedge-shaped structure makes it easy to observe the processing of the electrode. The thickness of the front part is 4mm and it holds cannulas 1mm to 3,5mm.

68U211

Fitting to the Stereotaxic Instruments with dual manipulators. It holds electrodes or cannulas of 2,1mm to 5,0mm in diameter. It can put them in distance of 0,3mm.



PAVAN SERVICE BVBA

Beyntellus 3 - B-2360 Oud-Turnhout - Tel. +32 (0)14 45 13 10 - admin@carfil.be - BTW BE 0420.681.872

www.carfil.be



CARFIL QUALITY
YOUR SUPPLIER FOR ALL LABORATORY NEEDS



68U217

It holds 3,5mm syringe pumps and cannular systems.



PAVAN SERVICE BVBA

Beyntellus 3 - B-2360 Oud-Turnhout - Tel. +32 (0)14 45 13 10 - admin@carfil.be - BTW BE 0420.681.872

www.carfil.be



CARFIL QUALITY
YOUR SUPPLIER FOR ALL LABORATORY NEEDS



68U208-M5

Suitable for dual cannula system of 5,0mm (OD). Working with manipulator on the stereotaxic, it implants the cannula into the targeted area of the brain.



PAVAN SERVICE BVBA

Beyntellus 3 - B-2360 Oud-Turnhout - Tel. +32 (0)14 45 13 10 - admin@carfil.be - BTW BE 0420.681.872

www.carfil.be



CARFIL QUALITY
YOUR SUPPLIER FOR ALL LABORATORY NEEDS



68U221 for Cannula with 0-3,5 mm size

Threaded locking mode is used to specifically hold the 3,5mm cannula system.



PAVAN SERVICE BVBA

Beyntellus 3 - B-2360 Oud-Turnhout - Tel. +32 (0)14 45 13 10 - admin@carfil.be - BTW BE 0420.681.872

www.carfil.be



CARFIL QUALITY
YOUR SUPPLIER FOR ALL LABORATORY NEEDS



68U209

This cannula holder employs the side clamp to secure the cannula or electrode more conveniently. It holds cannula with diameters from 2mm to 4,5mm.



PAVAN SERVICE BVBA

Beyntellus 3 - B-2360 Oud-Turnhout - Tel. +32 (0)14 45 13 10 - admin@carfil.be - BTW BE 0420.681.872

www.carfil.be