




CARFIL QUALITY
YOUR SUPPLIER FOR ALL LABORATORY NEEDS

Multi Face Mask



Features

- Allows to have 4 mice kept under anesthesia simultaneously
- Each mouse is placed with its snout directly into the face mask
- The multi face mask is delivered with stoppers that can seal the cones shut
- On the opposite of the openings for the mouse, the Multi facemask has four (4) connections for face masks with tubing. Thus making it possible to connect 4 rat facemasks with tubing. It is also possible to have 4 facemasks for mice connected if you want place the animals a bit further apart from eachother

Ordering number

Ordering number	Product	Details
180000136	Multi face mask	ID Supply: Ø 10 mm
		ID Exhaust: Ø 20 mm
		Overall dimensions: 264 x 115,5 x 33 mm
		Connection Supply hose: Ø 4 mm
		Connection Exhaust hose: Ø 22 mm
180000136R	Multi face Mask - extended for 4 rats	See Multi face mask plus 4 facemaks rat
180000136M	Multi face Mask - extended for 4 mice	See Multi face mask plus 4 facemask mouse



PAVAN SERVICE BVBA

Beyntellus 3 - B-2360 Oud-Turnhout - Tel. +32 (0)14 45 13 10 - admin@carfil.be - BTW BE 0420.681.872

www.carfil.be