

Face Mask for rat



Features

For inhalation anesthesia without respiration support, UNO has Face masks available

All models are to be connected to a supply system for anesthetic gas (like vaporizer and flowmeter or univentor). The second connection is to a gas exhaust system which safely removes the exhaled/ non used anesthetic gasses; like the UNO Gas Exhaust unit.

Face mask for rat

- With larger nose "cone" for gas supply
- To be used with rats (and mice depending on type of operation).
- Gas exhaust and cleanability as for the mouse face mask.
- * In combination with the UNO Gas Exhaust Unit.

Ordering number

Ordering number	Product	Details
180000237	Face Mask Rat	ID Supply: Ø 19 mm
		ID Exhaust: Ø 29 mm
		Overall dimensions: 81 x 33 x 33 mm
		Connection Supply hose: Ø 4 mm
		Connection Exhaust hose: Ø 22 mm





Fixator for face masks

The face masks are shaped with a flat bottom but sometimes you want to have it secured to a certain place. For this fixators are available.

- Made in Stainless Steel
- To keep the face mask in place, f.i. on a heating plate.