




**CARFIL QUALITY**  
YOUR SUPPLIER FOR ALL LABORATORY NEEDS

## Data Loggers DST Nano



### Eigenschappen

**DST Nano- worldcode2utf('39',0)s smallest implantable temperature logger.**

- Ideal for biomedical or animal welfare studies
- Constant and accurate temperature measurements
- Simple to use and cost effective
- Multiple animals in a cage (biomedical studies)
- Size only 17 mm x 6 mm, weighing only 1 g

The DST nano-T is the world's smallest temperature data logger. Recorded data is stored in the logger's internal memory with a real time clock reference for each measurement.

The DST nano-T is supported by the Mercury software and the Communication Box which serves as an interface between the logger and a PC. When logger is placed in the Communication Box the wireless communication between logger and box can be established. In the software, the user sets the start time, start date and sampling interval before starting the recorder. Up to seven different intervals can be set for the same measurement sequence. This is especially useful when more frequent measurements are needed at a certain time period. After recovering the DST, recorded data is uploaded in the supporting software where results are displayed both in graphic and tabular form. After retrieving the data, the DST can be re-programmed and reused as long as the battery lasts.

A set of Communication Box and Mercury software needs to be purchased with the first order.

### Examples of Application

DST nano-T was designed to be used as an implantable thermo logger for wild and captive animals. DST nano-T is made of biocompatible materials (alumina and biocompatible epoxy) and its size and shape are ideal for implantation in small animals. The logger can also be fitted in tight spaces, small packages or vials, for e.g. heat studies or temperature monitoring. With its waterproof housing and carefully chosen materials the DST nano-T can be used for collecting temperature data underwater or in other liquids (tolerates chemicals well).

The DST nano-T temperature logger is suitable for studies within:



**PAVAN SERVICE BVBA**

Beyntellus 3 - B-2360 Oud-Turnhout - Tel. +32 (0)14 45 13 10 - admin@carfil.be - BTW BE 0420.681.872

**www.carfil.be**



**CARFIL QUALITY**  
YOUR SUPPLIER FOR ALL LABORATORY NEEDS

- Laboratory animal science
- Preclinical research, stability and toxicology studies
- Quality control (e.g. packaging)
- Tagging of small fish and animals
- Any use where ultra small and lightweight temperature recorder is required

## Technische specificaties

|                                        |                                                                                           |
|----------------------------------------|-------------------------------------------------------------------------------------------|
| Sensors                                | Temperature                                                                               |
| Size                                   | 17 mm x Ø 6 mm                                                                            |
| Housing Material                       | Alumina and biocompatible epoxy                                                           |
| Weight in air/in water                 | 1 gram                                                                                    |
| Memory type                            | Non-volatile EEPROM                                                                       |
| Memory capacity                        | 43.476 measurements                                                                       |
| Memory capacity bytes                  | 65,534 Bytes / temp 1,5 byte                                                              |
| Memory management                      | Custom programming                                                                        |
| Data resolution                        | 12 bits                                                                                   |
| Temperature resolution                 | 0,032°C                                                                                   |
| Temperature accuracy                   | +/- 0,2°C                                                                                 |
| Temperature range                      | 5 to 45°C                                                                                 |
| Temperature response time              | Time constant (63%) reached in 5sec                                                       |
| Data retention                         | 10 years                                                                                  |
| Clock                                  | Real time clock, accuracy +/- 1 min/month                                                 |
| Sampling Interval                      | User specified in seconds, minutes or hours                                               |
| Number of different sampling intervals | 1 to 7                                                                                    |
| Communications                         | Communication box, RS232C 9 pin serial or USB                                             |
| Battery life                           | 14 months - with a sampling interval of 10 min at room temp.<br>Non-replaceable batteries |

## Bestelnummer

180000093

DST nano-T temperature recorder



**PAVAN SERVICE BVBA**

Beyntellus 3 - B-2360 Oud-Turnhout - Tel. +32 (0)14 45 13 10 - admin@carfil.be - BTW BE 0420.681.872

**www.carfil.be**



**CARFIL QUALITY**  
YOUR SUPPLIER FOR ALL LABORATORY NEEDS



**PAVAN SERVICE BVBA**

Beyntellus 3 - B-2360 Oud-Turnhout - Tel. +32 (0)14 45 13 10 - [admin@carfil.be](mailto:admin@carfil.be) - BTW BE 0420.681.872

**[www.carfil.be](http://www.carfil.be)**