



# **Tether Infusion Sets & Extension Sets**



# **TIS Series**

### **Blunt Needles Tether Infusion Sets**

Tether Infusion Sets feature a blunt needle bonded to a Tygon extension and are used to connect a catheter to the infusion source through a tether. They can be ordered in a variety of needle gauges with custom needle and extension tubing lengths and are fitted with a female luer unless otherwise specified. Tether Infusion Sets ship sterile (EtO) in boxes of 20.

### **Tish Series**

### **Huber Needle Tether Infusion Sets**

Huber Tether Infusion Sets feature a Huber needle bonded to a Tygon extension and are used to access a Vascular Access Port. The infusion line passes from the V-A-P through a tether system to the infusion source. They can be ordered in a variety of needle gauges with custom needle and extension tubing lengths and are fitted with a female luer unless otherwise specified. Huber Tether Infusion Sets ship sterile (EtO) in boxes of 20.

# Ordering Information

Needle Gauge	Tubing Length	Needle Length							
		1/2" (12,7 mm)	5/8" (15,88 mm)	3/4" (19,05 mm)	1" (25,4 mm)				
TISH Series - Huber Tether Infusion Sets - with a Huber needle, female luer & clamp									
22	45 cm	TISH22-5-18	TISH22-625-18	TISH22-75-18	TISH22-100-18				
Beyntellus 3 -		<sup></sup> el. +32 (0)14 45 13 10 - a	dmin@carfil.be - BTW BE	0420.681.872	www.carfil.b				



### **CARFIL QUALITY** YOUR SUPPLIER FOR ALL LABORATORY NEEDS

22	90 cm	TISH22-5-36	TISH22-625-36	TISH22-75-36	TISH22-100-36				
20	45 cm	TISH20-5-18	TISH20-625-18	TISH20-75-18	TISH20-100-18				
20	90 cm	TISH20-5-36	TISH20-625-36	TISH20-75-36	TISH20-100-36				
TIS Serie - Tether Infusion Sets - with a Blunt needle, female luer & clamp									
22	45 cm	TIS22-5-18	TIS22-625-18	TIS22-75-18	TIS22-100-18				
22	90 cm	TIS22-5-36	TIS22-625-36	TIS22-75-36	TIS22-100-36				
20	45 cm	TIS20-5-18	TIS20-625-18	TIS20-75-18	TIS20-100-18				
20	90 cm	TIS20-5-36	TIS20-625-36	TIS20-75-36	TIS20-100-36				

# ES Series Extensions Sets

Extensions Sets are available in a variety of materials, including Polyurethane, Tygon and PE-PVC Coextruded tu- bing in custom lengths and configurations. They ship sterile (EtO) in boxes of 10 with male and female luers and a clamp unless otherwise specified. When ordering Extension Sets, specify the tubing material, length and luer configuration. Examples: ES-60M/M is a Tygon extension set 60" long with a clamp & male luers on both ends.

ES-7P-48M/F is a Polyurethane extension set 48" long with male and female luers.

### Luered Extension Lines - Catheter to Swivel connection

Polyurethane or Tygon Extension Lines with catheter connectors or male luers fit a wide range of swivels. These extension sets connect the catheter from the animal to a swivel. The connection is made with a stainless steel catheter connector or a luer on the end of the catheter.

### Luered Extension Lines - Male to Female

Polyurethane or Tygon Extension Lines with various luer configurations and lengths: male-male, male-female or female-female.

### Luered Extension Lines - Swivel to Pump connection

Polyurethane or Tygon Extension Lines in a variety of diameters and with female luers to fit most swivel gauges





www.carfil.be